



## DRAWING CONTENTS

**C1 COVER SHEET**  
**GN1 GENERAL NOTES**  
**GN2 GENERAL NOTES**

**SECTION A**  
**A1 EQUIPMENT LAYOUT**  
**A2 EQUIPMENT ELEVATIONS**

**SECTION S**  
**S1 STRUCTURAL LAYOUT**

**SECTION E**  
**E1 ELECTRICAL LAYOUT**  
**E2 ELECTRICAL SCHEMATIC**  
**E3 ELECTRICAL DETAILS**  
**E4 CABLE LENGTH DIAGRAM (FOR INTERNAL USE ONLY)**

HARRISBURG HOSPITAL		REV	DATE	REVISED SHEET(S)	INT
(AQUILION – RXL)					
111 S FRONT ST HARRISBURG, PA 17101					
THESE TOSHIBA PLANS ARE FOR INFORMATIONAL PURPOSES ONLY AND SHALL NOT BE USED FOR ANY PURPOSE OTHER THAN THAT AGREED UPON BETWEEN TOSHIBA AND THE CUSTOMER. THESE SITE PLANS ARE NOT TO BE USED FOR CONSTRUCTION PURPOSES.					
DATE:		06-19-13			
SCALE:		NOT TO SCALE			
PLANNER:		C.B.S.			
QUOTE:		N/A			
PROJECT NO.		<b>130013698CTF</b>			
		<b>C1</b>			

## GENERAL NOTES

## ELECTRICAL REQUIREMENTS FOR AQUILION

**CUSTOMER / CONTRACTOR IS RESPONSIBLE FOR THE FOLLOWING UNLESS OTHERWISE NOTED.**

**GENERAL**

A. TOSHIBA RESERVES THE RIGHT TO CHANGE THESE DESIGNS AND SPECIFICATIONS WITHOUT NOTICE.

B. THE CUSTOMER/CONTRACTOR IS RESPONSIBLE TO ENSURE THAT ALL APPLICABLE FEDERAL, STATE, AND LOCAL CODES AND ORDINANCES ARE COMPLIED WITH.

C. PRIOR TO EQUIPMENT DELIVERY AND INSTALLATION, THE SITE MUST BE 100% COMPLETE, CLEAN AND FREE OF DUST. CUSTOMER/CONTRACTOR AND TOSHIBA INSTALLATION PROJECT MANAGER MUST COMPLETE A SITE WALK THROUGH 1 WEEK PRIOR TO DELIVERY AND DETERMINE ACCEPTABILITY FOR DELIVERY.

D. ANY CABINETY THAT MAY BE REQUIRED TO HOUSE VIDEO RECORDERS, MONITORS KEYBOARDS, OR OTHER ANCILLARY EQUIPMENT SHALL BE SUPPLIED AND INSTALLED BY CUSTOMER/CONTRACTOR.

E. PROVIDE ADEQUATE VENTILATION WITHIN CABINETY AND INSTALL AXIAL FANS ON THE TOP, SIDE, OR BACK OF CABINETS, IF REQUIRED.

F. THESE TOSHIBA SITE PLANS DO NOT INDICATE EQUIPMENT REQUIREMENTS FOR ITEMS NOT SOLD BY TOSHIBA SUCH AS, PHYSIOLOGICAL MONITORS, LASER CAMERAS, INJECTORS, ETC. SPECIFICATIONS FOR THOSE ITEMS MUST BE OBTAINED FROM THE VENDOR AND INCLUDED IN THE DESIGN TOTALS.

G. DESIGN, FABRICATE, AND INSTALL MEDICAL GAS PEDESTAL, IF REQUIRED. CONSULT WITH TOSHIBA INSTALLATION PROJECT MANAGER FOR SUITABLE LOCATIONS.

H. CUSTOMER/CONTRACTOR SHALL BE RESPONSIBLE FOR PROVIDING AN OPERATING PHONE IN THE CONTROL ROOM AT THE TIME TOSHIBA EQUIPMENT INSTALLATION BEGINS.

I. CUSTOMER/CONTRACTOR SHALL BE RESPONSIBLE FOR PROVIDING ADEQUATE LIGHTING FOR SERVICING OF EQUIPMENT IN ALL AREAS OF THE INSTALLATION.

J. THE CUSTOMER/CONTRACTOR SHALL BE RESPONSIBLE FOR ANY AND ALL COSTS REQUIRED FOR THE ENGINEERING AND/OR REMOVAL OF ANY HAZARDOUS MATERIALS SUCH AS ASBESTOS.

K. CUSTOMER/CONTRACTOR SHALL SUPPLY AND INSTALL MATERIALS AND OTHER FEATURES SPECIFIED IN THE TOSHIBA SITE PLANS. CUSTOMER/CONTRACTOR SHALL SUPPLY AND INSTALL ALL COUNTERTOPS, SINKS, CASE WORK AND CABINETS SPECIFIED IN THE TOSHIBA SITE PLANS.

**PLUMBING**

L. PLUMBING IS NOT REQUIRED FOR THIS TOSHIBA EQUIPMENT.

M. IT IS RECOMMENDED THAT A SINK BE PROVIDED FOR USE BY PERSONNEL.

**SITE CONDITIONS**

N. DIMENSIONS TO WALLS AND OR OTHER ROOM FEATURES, EXCEPT FOR NOTED COLUMN AND BEAM CENTER LINES SHALL BE FROM FINISHED SURFACES.

O. CT GANTRY SHOULD NOT BE INSTALLED WITHIN 0.5 MAGNETIC GAUSS FIELD.

P. THE WINDOW FOR MONITORING THE SCAN ROOM SHOULD BE IN FRONT OF OR ON THE SIDE OF THE CONSOLE DESK. THE LOWEST WINDOW FRAME SHOULD BE 36" ABOVE THE FLOOR FOR EASY PATIENT MONITORING.

Q. A DOOR BETWEEN THE SCAN AND CONTROL ROOM IS RECOMMENDED.

R. THE INSTALLATION ALTITUDE SHOULD BE NO MORE THAN 3,280 FT. (1,000 M) ABOVE SEA LEVEL. PRIOR CONSULTATION IS REQUIRED FOR INSTALLATIONS HIGHER THAN 3,280 FT. (1,000 M).

**NETWORKING REQUIREMENTS**

S. NETWORK REQUIREMENTS WILL VARY BY SITE. TOSHIBA REPRESENTATIVE WILL REQUIRE DICOM DEVICE INFORMATION. ADDITIONAL I.P. ADDRESSES, AND I.T. DEPARTMENT CONTACT INFORMATION PRIOR TO INSTALLATION.

**TRANSPORT REQUIREMENTS**

T. EQUIPMENT INGRESS ROUTE MUST BE CHECKED PRIOR TO EQUIPMENT DELIVERY TO ENSURE THE LARGEST AND HEAVIEST ITEMS OF EQUIPMENT CAN BE ACCOMMODATED, PRIOR TO EQUIPMENT DELIVERY. DIMENSIONS OF DOORWAYS SHOULD BE NO LESS THAN 4'-0" CLEAR IN WIDTH. CONTACT THE TOSHIBA INSTALLATION PROJECT MANAGER FOR DETAILS PERTAINING TO THE LARGEST AND HEAVIEST COMPONENTS FOR THIS INSTALLATION (SEE DETAIL 4, SHEET GN2).

02-13-13

02-13-13

02-13-13

02-13-13

04-09-13 DATE: 06-19-13

04-09-13

**FOR REFERENCE ONLY. NOT TO BE USED FOR CONSTRUCTION PURPOSES.**

DATE: 06-19-13

SCALE: NOT TO SCALE

PLANNER: C.B.S.

QUOTE: N/A


PROJECT NO.  
**130013698CTF**

# GN1

CUSTOMER TO PROVIDE THE NECESSARY HVAC REQUIREMENTS FOR THE TOSHIBA EQUIPMENT TO OPERATE PROPERLY.

AMBIENT TEMPERATURE SHOULD BE 68°–74° F  
WITH EQUIPMENT HEAT LOADS (SEE EQUIPMENT LEGEND SHEET A1)  
HUMIDITY RANGE OF 40–70% NON–CONDENSING

- A. STATED AMBIENT TEMPERATURE IS TO BE PROVIDED AND MAINTAINED AS SPECIFIED. ALL CALCULATIONS ARE TO UTILIZE TOSHIBA PROVIDED HEAT OUTPUT SPECIFICATIONS OF EQUIPMENT.
- B. A MINIMUM OF 10 AIR CHANGES PER HOUR IS SUGGESTED, CONSULT LOCAL CODE.
- C. AIR SUPPLY DUCTS SHOULD NOT BE PLACED DIRECTLY OVER EXAMINATION TABLES FOR PATIENT COMFORT (SEE DETAIL 2, SHEET GN2).
- D. EQUIPMENT IN ENCLOSED SPACES SUCH AS EQUIPMENT ROOMS, TRANSFORMER CLOSETS AND COMPUTER ROOMS MUST BE PROVIDED WITH ADEQUATE VENTILATION.
- E. THE AIRFLOW THROUGH TOSHIBA EQUIPMENT CABINETS IS FROM BOTTOM TO TOP.
- F. WHERE POSSIBLE, AIR CONDITIONING SUPPLY OUTLETS SHOULD BE LOCATED AT FLOOR LEVEL. NO AIR CONDITIONING OUTLET SHOULD BE WITHIN THE EXCLUSION ZONE SHOWN BELOW AND AT NO TIME SHOULD THE CT SYSTEM BE EXPOSED TO DIRECT AIRFLOW.
- G. RETURN GRILLES ARE TO BE INSTALLED IN THE CEILING.
- H. A/C SUPPLY OUTLET TO BE PROVIDED BY CUSTOMER AT FLOOR LEVEL AT CONTROL ROOM DESK.
- I. DUE TO HEAT GENERATED BY THE "CPU" UNIT, ADDITIONAL VENTILATION IN THE CONTROL AREA IS REQUIRED. CUSTOMER/CONTRACTOR PROVIDED FAN(S) MAY BE NECESSARY BELOW THE DESKTOP FOR TECHNICIAN COMFORT. THE "CPU" UNIT SHOULD NOT BE ENCLOSED IN CASEWORK.



- J. IN GENERAL, THE SCANNING ROOM MUST BE PROVIDED WITH AN INDEPENDENT AIR CONDITIONING SYSTEM. EVEN IF THE ROOM IS MAINTAINED WITHIN THE PERMISSIBLE TEMPERATURE RANGE, GRADUAL TEMPERATURE SHIFTS (FOR EXAMPLE, A SLOW INCREASE IN ROOM TEMPERATURE FROM MORNING TO EVENING) MAY ADVERSELY AFFECT SYSTEM PERFORMANCE. THEREFORE, THE ROOM TEMPERATURE MUST BE KEPT UNDER CONSTANT CONTROL (WITHIN ±3°F) AS SHOWN IN THE ABOVE FIGURE.
- K. THE AIR CONDITIONING SYSTEM IN THE SCANNING ROOM MUST BE INSTALLED SO THAT THE CT SYSTEM IS NOT EXPOSED TO DIRECT AIRFLOW. FAILURE TO DO SO MAY CAUSE THE TEMPERATURE INSIDE THE CT SYSTEM TO FLUCTUATE, POSSIBLY AFFECTING THE DISPLAYED IMAGES ADVERSELY.

1 HVAC REQUIREMENTS


SCALE: NOT TO SCALE

02–13–13

2 SUPPLY OUTLET EXCLUSION ZONE

SCALE: 1/4" = 1'–0"

02–13–13




- A. CEILING MOUNTED EQUIPMENT MUST BE POSITIONED TO AVOID INTERFERENCE WITH GANTRY.
- B. EQUIPMENT IS TO BE A MINIMUM OF 10" ABOVE THE RAISED GANTRY COVER (SEE DETAIL 1, SHEET A2).
- C. OVERHEAD COUNTERPOISE SYSTEMS CAN BE INSTALLED OVER GANTRY ISOCENTER IF THE PLATE IS MOUNTED AT A HEIGHT SUCH THAT THE BOTTOM OF THE POST DOES NOT INTERFERE WITH THE GANTRY (HEIGHT OF GANTRY + 10" CLEARANCE + POST LENGTH = MOUNTING PLATE HEIGHT ABOVE FINISHED FLOOR).


3 CEILING MOUNTED EQUIPMENT

SCALE: 1/4" = 1'–0"

02–13–13



PRIOR TO GANTRY DELIVERY, CHECK THE DELIVERY ENTRANCE, CORRIDOR WIDTH, ELEVATOR CAPACITY, ETC. TO ENSURE EASE OF DELIVERY. COMPENSATION FOR ADDITIONAL REQUIREMENTS OF RIGGING EQUIPMENT USED (I.E. SIZE AND MANEUVERABILITY) MUST BE CONSIDERED WHEN REVIEWING THE FOLLOWING.



MINIMUM DELIVERY REQUIREMENTS

MINIMUM HEIGHT OF ENTRY WAY DOORS: 6'–9"

MINIMUM WIDTH OF ENTRY WAY DOORS: 3'–7 1/2"

MINIMUM WIDTH OF CORRIDOR: 6'–7"

4 RECOMMENDED GANTRY DELIVERY REQUIREMENTS

SCALE: NOT TO SCALE

02–13–13

HARRISBURG  
HOSPITAL

(AQUILION – RXL)

111 S FRONT ST  
HARRISBURG, PA 17101

THESE TOSHIBA PLANS ARE FOR INFORMATIONAL PURPOSES ONLY AND SHALL NOT BE USED FOR ANY PURPOSE OTHER THAN THAT AGREED UPON BETWEEN TOSHIBA AND THE CUSTOMER. THESE SITE PLANS ARE NOT TO BE USED FOR CONSTRUCTION PURPOSES.

DATE: 06–19–13

SCALE: NOT TO SCALE

PLANNER: C.B.S.

QUOTE: N/A


PROJECT NO.  
130013698CTF

GN2

TOSHIBA  
Leading Innovation >>>

FOR REFERENCE ONLY. NOT TO BE USED FOR CONSTRUCTION PURPOSES.





## SCATTER RADIATION

UNITS: MICROGRAY (PER 100 mAs)

THE AMOUNT OF SCATTERED RADIATION DURING SCANNING IS SHOWN IN THE EQUIPMENT LAYOUT. REFER TO THIS DATA WHEN SCANNING TO MINIMIZE X-RAY EXPOSURE. NOTE THAT THE AMOUNT OF SCATTERED RADIATION DIFFERS DEPENDING ON THE PATIENT; THEREFORE, THE DATA SHOWN SHOULD BE USED ONLY AS A GUIDE.


120kV / 100mA / 1.0s / M / 8mm x 4 / 320mm diameter PMMA phantom

THE ABOVE CALCULATION IS BASED ON A  
MULTI-SLICE EXPOSURE.

THE CUSTOMER/CONTRACTOR IS RESPONSIBLE  
FOR HAVING SHIELDING CALCULATIONS  
PREPARED BY A LICENSED RADIATION PHYSICIST

TOSHIBA REPRESENTATIVE WILL PROVIDE THE  
VERTICAL PLANE SCATTER DIAGRAM UPON  
REQUEST.

REVISÉ: 02-13-13



# SITE PLAN APPROVAL


IN ORDER TO USE THIS SET OF FINAL SITE PLANS, A CUSTOMER SIGNATURE IS REQUIRED BELOW. THE CUSTOMER'S SIGNATURE DEMONSTRATES ACCEPTANCE OF THE LAYOUT SHOWN AND ALL STATED SPECIFICATIONS.

CUSTOMER:	DATE:
SALES:	DATE:
I.P.M.:	DATE:

<div style="text-align: right;"> <b>TOSHIBA</b>  <b>Leading Innovation &gt;&gt;&gt;</b> </div>					
REV	DATE	DESCRIPTION	INT		
<b>HARRISBURG HOSPITAL</b>  (AQUILION – RXL)  111 S FRONT ST HARRISBURG, PA 17101		<p>THESE TOSHIBA PLANS ARE FOR INFORMATIONAL PURPOSES ONLY AND SHALL NOT BE USED FOR ANY PURPOSE OTHER THAN THAT AGREED UPON BETWEEN TOSHIBA AND THE CUSTOMER. THESE SITE PLANS ARE NOT TO BE USED FOR CONSTRUCTION PURPOSES.</p>			
DATE:		06-19-13			
SCALE:		1/4" = 1'-0"			
PLANNER:		C.B.S.			
QUOTE:		N/A			
PROJECT NO. <b>130013698CTF</b>					
<b>A1</b>					

**FOR REFERENCE ONLY. NOT TO BE USED FOR CONSTRUCTION PURPOSES.**






STRUCTURAL LAYOUT



NOTE:  
ANY FLOOR COVERING SHOULD BE REMOVED WITHIN THE AREA SHOWN IN DETAIL 1 EXPOSING CONCRETE PRIOR TO POURING THE EPOXY PAD. THE AREA FROM WHICH THE FLOOR COVERING HAS BEEN REMOVED MUST BE CLEAN AND FREE OF RESIDUE.

STRUCTURAL LEGEND		
ITEM	ITEM DESCRIPTION SUPPLIED AND INSTALLED BY TOSHIBA	REF.
A	BASE LAYOUT FOR AQUILON	1 S1



DETAIL NOTES

AN EPOXY PAD MUST BE POURED TO PROVIDE A LEVEL SURFACE. AN EPOXY PAD KIT WILL BE PROVIDED WITH THE PRE-INSTALLATION MATERIALS. CONTACT YOUR INSTALLATION PROJECT MANAGER FOR ADDITIONAL INFORMATION.


APPLY EPOXY RESIN TO THE AREA INDICATED ABOVE.

- EPOXY RESIN: SPECIFIC GRAVITY 1.2  
HARDENER: SPECIFIC GRAVITY 1.2  
MIXTURE RATIO:PER MANUFACTURER'S SPECIFICATIONS  
CURING TIME: 36-48 HOURS (AT AN AMBIENT TEMPERATURE OF APPROXIMATELY 75°F)
- SINCE THE EPOXY RESIN TAKES 36-48 HOURS (DEPENDING ON AMBIENT TEMPERATURE) TO CURE, THIS WORK MUST BE COMPLETED BEFORE THE SYSTEM IS TO BE CARRIED IN.
- THE ACCURACY OF THE DIMENSIONS INDICATED IN THE ABOVE FIGURE MUST BE WITHIN THE RANGE OF  $\pm 1/4$ " AS MEASURED WITH A TAPE MEASURE.

1 BASE EPOXY DETAIL

SCALE: 1/2" = 1'-0"

02-13-13



2 BASE ANCHOR BOLT DETAIL

SCALE: 1/2" = 1'-0"

02-13-13

TOSHIBA

Leading Innovation >>>

INT	DESCRIPTION	DATE	REV

HARRISBURG HOSPITAL

(AQUILON - RXL)

111 S FRONT ST  
HARRISBURG, PA 17101

THESE TOSHIBA PLANS ARE FOR INFORMATIONAL PURPOSES ONLY AND SHALL NOT BE USED FOR ANY PURPOSE OTHER THAN THAT AGREED UPON BETWEEN TOSHIBA AND THE CUSTOMER. THESE SITE PLANS ARE NOT TO BE USED FOR CONSTRUCTION PURPOSES.

DATE: 06-19-13

SCALE: 1/4" = 1'-0"


PLANNER: C.B.S.

QUOTE: N/A

PROJECT NO.  
**130013698CTF**


S1


FOR REFERENCE ONLY. NOT TO BE USED FOR CONSTRUCTION PURPOSES.




**NOTE:**  
GROMMETED OPENINGS ARE SHOWN FOR  
REFERENCE ONLY. VERIFY SIZE AND LOCATION  
WITH TOSHIBA REPRESENTATIVE.

ELECTRICAL DUCT LEGEND		
ITEM	ITEM DESCRIPTION SUPPLIED AND INSTALLED BY CUSTOMER / CONTRACTOR	REF.
WD1	EXISTING SURFACE MOUNTED WALL DUCT, W/(3) EQUALLY PARTITIONED COMPARTMENTS THROUGHOUT & REMOVABLE ACCESS COVERS. MOUNTED 3" OR A.F.F. TO BOTTOM OF DUCT.	
WD2	EXISTING SURFACE MOUNTED METAL WALL DUCT W/(3) EQUAL PARTITIONED COMPARTMENTS THROUGHOUT & REMOVABLE ACCESS COVERS. MOUNTED 3" A.F.F. TO BOTTOM OF DUCT. OPEN TO "WD1" VIA (3) 3" CHASE NIPPLES.	
VR1	EXISTING SURFACE MOUNTED RISER DUCT, W/(3) EQUALLY PARTITIONED COMPARTMENTS THROUGHOUT & REMOVABLE ACCESS COVERS. FROM "PD" TO "WD2".	
VR2	EXISTING SURFACE MOUNTED RISER DUCT, W/(3) EQUALLY PARTITIONED COMPARTMENTS THROUGHOUT & REMOVABLE ACCESS COVERS. FROM "CD" TO "WD1".	
CD	EXISTING DUCT MOUNTED ABOVE FINISHED CEILING, W/(3) EQUALLY PARTITIONED COMPARTMENTS THROUGHOUT & REMOVABLE ACCESS COVERS (THICKNESS SHOULD BE ABLE TO SUPPORT MIN. OF 200 LBS.) ARE REQUIRED.	


<div style="text-align: center;">  <h1>TOSHIBA</h1> <p>Leading Innovation &gt;&gt;&gt;</p> </div>						INT
<b>HARRISBURG HOSPITAL</b>  (AQUILION — RXL)  111 S FRONT ST HARRISBURG, PA 17101						
THESE TOSHIBA PLANS ARE FOR INFORMATIONAL PURPOSES ONLY AND SHALL NOT BE USED FOR ANY PURPOSE OTHER THAN THAT AGREED UPON BETWEEN TOSHIBA AND THE CUSTOMER; THESE SITE PLANS ARE NOT TO BE USED FOR CONSTRUCTION PURPOSES.						
DATE:		06-19-13				
SCALE:		$\frac{1}{4}" = 1'-0"$				
PLANNER:		C.B.S.				
QUOTE:		N/A				
PROJECT NO. <b>130013698CTF</b>						
<b>E1</b>						



NOTE:  
REUSE ANY EXISTING DUCT, J-BOXES OR  
CONDUITS IF APPLICABLE.



ELECTRICAL SCHEMATIC  
(PROVIDED FOR REFERENCE PURPOSES ONLY)



CABLE KEY	
	IN/UNDER FLOOR
	OVER CEILING
	CONTRACTOR DETERMINED

CONDUIT SCHEDULE

CONTRACTOR CONDUIT REFERENCE						CABLE REFERENCE			
RUN NO.	CONDUIT (POINT TO POINT)		CONDUIT (ROUTING)	CONDUIT (DIAMETER)	CONDUIT (MAX LENGTH)	CABLE (POINT TO POINT)		CABLE LENGTH (USABLE)	CABLES (SUPPLIED BY)
①	WD1	GANT	UNDER FLOOR	(2) 2 1/2"	20'-0"	PD	GANT	SEE RUN "B" DETAIL (1/E4)	TOSHIBA
②	GANT	WD1	UNDER FLOOR	2 1/2"	20'-0"	GANT	CPU	SEE RUN "D" DETAIL (1/E4)	TOSHIBA
③	GANT	WD1	UNDER FLOOR	(2) 3"	20'-0"	GANT	CPU	SEE RUN "E" DETAIL (1/E4)	TOSHIBA
④	SPK2	CD	OVER CEILING	1/2"	25'-0"	SPK2	CPU	SEE RUN "G" DETAIL (1/E4)	TOSHIBA
⑤	WD1	GANT	UNDER FLOOR	2 1/2"	20'-0"	REC	GANT	SEE RUN "L" DETAIL (1/E4)	TOSHIBA
⑥	WL	CD	OVER CEILING	PER MANUFACTURER	PER MANUFACTURER	WL	GANT	PER MANUFACTURER	CONTRACTOR
⑦	GANT	WD1	UNDER FLOOR	1"	20'-0"	GANT	INV	SEE RUN "H" DETAIL (1/E4)	TOSHIBA

NOTE:  
A. CONDUITS SUPPLIED/INSTALLED BY CUSTOMER/CONTRACTOR.  
B. ALL CONDUIT RUNS MUST TAKE THE SHORTEST MOST DIRECT ROUTE POSSIBLE.  
C. \* IF RUN IS GREATER THAN LENGTH SHOWN, CUSTOMER/CONTRACTOR TO PROVIDE CABLES.  
REFER TO DETAIL 1, SHEET E3.  
D. CONDUIT IS NOT TO BE RUN IN SUCH A MANNER THAT WILL EXCEED CONDUIT MAXIMUM LENGTH AS SHOWN IN THE SCHEDULES.

TOSHIBA  
Leading Innovation >>>

INT

DESCRIPTION

DATE

REV

HARRISBURG  
HOSPITAL

(AQUILION – RXL)

111 S FRONT ST  
HARRISBURG, PA 17101

THESE TOSHIBA PLANS ARE FOR INFORMATIONAL PURPOSES ONLY AND SHALL NOT BE USED FOR ANY PURPOSE OTHER THAN THAT AGREED UPON BETWEEN TOSHIBA AND THE CUSTOMER. THESE SITE PLANS ARE NOT TO BE USED FOR CONSTRUCTION PURPOSES.

DATE: 06-19-13

SCALE: 1/4" = 1'-0"

PLANNER: C.B.S.

QUOTE: N/A




PROJECT NO.  
130013698CTF

E2

FOR REFERENCE ONLY. NOT TO BE USED FOR CONSTRUCTION PURPOSES.




POWER QUALITY REQUIREMENTS AQILION		
SUPPLY CONFIGURATION:	3 PHASE, 3 WIRE POWER, AND GROUND, DELTA OR WYE (SEE NOTE A)	
NOMINAL LINE VOLTAGE:	480 VAC, 60 HZ (SEE NOTE B)	
LINE VOLTAGE VARIATION:	±10% STEADY-STATE INCLUDING SAGS AND SURGES	
LINE VOLTAGE BALANCE:	2% MAXIMUM OF NOMINAL VOLTAGE BETWEEN PHASES.	
FREQUENCY VARIATION:	±1 HZ	
HARMONIC DISTORTION:	3% STEADY-STATE, 5% FOR SHORT PERIODS (1 MINUTE OR LESS)	
GROUND CONDUCTOR IMPEDANCE:	0.1 OHMS @ 60 HZ, TO NEUTRAL-GROUND BONDING POINT (SEE NOTE D)	
STANDARD TRANSFORMER CAPACITY:	150 KVA	
MAXIMUM SYSTEM DEMAND:	100 KVA (IMAGING)	
CONDUCTOR SIZES (SEE NOTE E) FOR 1.5% IMPEDANCE OF BRANCH CONDUCTORS (20°C)		
CONDUCTOR SIZE	LENGTH	BREAKER FRAME SIZE
#2 AWG	218 FT.	200 A
#1 AWG	303 FT.	200 A
1/0 AWG	377 FT.	200 A
2/0 AWG	444 FT.	200 A
CIRCUIT BREAKER SIZE: (SEE NOTE F)	100 A	
MOMENTARY MAXIMUM CURRENT:	150 A	
MAXIMUM VOLTAGE DROP:	24.0 V	
% REGULATION:	5.0%	



## STANDARD POWER QUALITY NOTES

- A GROUNDED NEUTRAL POWER SOURCE IS REQUIRED TO ASSURE RELIABLE EQUIPMENT OPERATION. THE NEUTRAL CONDUCTOR MAY NOT BE USED FOR A PARTICULAR SYSTEM.
- IN CASES WHERE MULTIPLE VOLTAGES ARE PERMITTED, THE PREFERRED SYSTEM VOLTAGE IS SPECIFIED.
- DUE TO THE HIGH INSTANTANEOUS POWER OF MEDICAL IMAGING SYSTEMS, USE THE HIGHEST AVAILABLE VOLTAGE SOURCE. ENSURE THAT LOWER VOLTAGE SOURCES ARE DERIVED DIRECTLY FROM THE SERVICE ENTRANCE OF THE FACILITY.
- GROUND CONDUCTORS ARE REQUIRED TO BE THE SAME SIZE AS THE PHASE CONDUCTORS UNLESS A LARGER SIZE IS REQUIRED BY CODE.
- ALL FEEDER AND BRANCH CIRCUIT CONDUCTORS MUST BE COPPER - ALUMINUM IS NOT PERMITTED.
- IF THE EQUIPMENT CIRCUIT BREAKER IS NOT LOCATED IN THE CONTROL AREA, A SHUNT TRIP BREAKER MUST BE USED IN ORDER TO COMPLY WITH N.E.C. 517-72(B). A PUSH-BUTTON TO OPERATE THE SHUNT TRIP MUST BE LOCATED IN THE CONTROL AREA.
- A SEPARATE CIRCUIT, FED FROM THE FACILITY RADIOLOGY PANEL OR A MAIN SERVICE PANEL IS REQUIRED. USE OF A SUB PANEL WITH LOADS SUCH AS ELEVATORS, HVAC, MOTORS, ETC., IS NOT PERMITTED.
- DEVICES SUCH AS UNINTERRUPTIBLE POWER SUPPLIES, POWER CONDITIONERS, VOLTAGE REGULATORS, AND FILTERS MAY BE INCOMPATIBLE WITH THIS IMAGING EQUIPMENT. CONSULT YOUR TOSHIBA SERVICE REPRESENTATIVE PRIOR TO PURCHASING OR INSTALLING THESE DEVICES.
- THE MAINS POWER GROUND CONDUCTOR IS TO BE RUN WITH THE POWER PHASE CONDUCTORS. THE GROUNDS TO BUILDING STEEL OR EARTH GROUND ARE NOT TO BE RUN WITH THE PHASE CONDUCTORS.



## 2 "CB" / "PD" WIRING DETAIL

SCALE: NOT TO SCALE

02-25-13

## 3 WARNING LIGHT DETAIL


SCALE: NOT TO SCALE

04-09-13

## 4 TYPICAL "PD" J-BOX MOUNTING

SCALE: 1/2" = 1'-0"

02-13-13



## 5 TYPICAL DUCT DETAIL WITH WALL DUCT / J-BOX / VERTICAL RISER

SCALE: NOT TO SCALE

04-09-13

## 1 POWER REQUIREMENTS

SCALE: NOT TO SCALE

03-06-13

REV	DATE	DESCRIPTION	INT

**HARRISBURG  
HOSPITAL**

(AQILION - RXL)

111 S FRONT ST  
HARRISBURG, PA 17101

THESE TOSHIBA PLANS ARE FOR INFORMATIONAL PURPOSES ONLY AND SHALL NOT BE USED FOR ANY PURPOSE OTHER THAN THAT AGREED UPON BETWEEN TOSHIBA AND THE CUSTOMER. THESE SITE PLANS ARE NOT TO BE USED FOR CONSTRUCTION PURPOSES.


DATE: 06-19-13

SCALE: NOT TO SCALE

PLANNER: C.B.S.

QUOTE: N/A

PROJECT NO.  
**130013698CTF**



NOTE:  
A. THE FOLLOWING CABLE DIAGRAM IS FOR REFERENCE ONLY AND MAY NOT CONTAIN THE EXACT EQUIPMENT SHOWN IN THE ATTACHED DRAWING PACKAGE.  
B. INSTALLATION PROJECT MANAGER IS RESPONSIBLE FOR ENSURING THAT CABLE LENGTHS MEET THE SITE CONDITIONS. THE INSTALLATION PROJECT MANAGER IS RESPONSIBLE FOR ORDERING THE REQUIRED CABLES TO MEET EXISTING/PROPOSED SITE CONDITIONS.  
C. ALL CABLE CONNECTIONS SHOWN ARE MAXIMUM LENGTH AND CANNOT BE EXTENDED.  
D. THIRD PARTY ITEM CONNECTIONS TO BE VERIFIED WITH INSTALLATION PROJECT MANAGER.  
E. POWER AND SIGNAL CABLES MUST BE RUN IN SEPARATE CONDUIT. REFER TO SHEET E2.

RUN LTR.	CONNECTION (POINT TO POINT)	CABLE LENGTH (MAX. USABLE)
A	POWER SYSTEM PD	74'-9"
		74'-9"
B	PD GANT	75'-1"
		75'-1"
		74'-9"
C	PD CPU	77'-5"
		75'-1"
D	GANT CPU	73'-9"
		74'-9"
E	GANT CPU	73'-9"
F	CPU SPK1	28'-6"
G	SPK2 CPU	48'-2"
H	GANT INV	75'-0"
I	CPU INV	75'-0"
J	REC CPU	62'-0"
K	REC CPU	63'-7"
		63'-7"
		63'-7"
L	GANT REC	73'-9"

TOSHIBA

Leading Innovation >>>

INT											
DESCRIPTION											
DATE											
REV											
HARRISBURG HOSPITAL		(AQUILION – RXL)	111 S FRONT ST	HARRISBURG, PA 17101							
THESE TOSHIBA PLANS ARE FOR INFORMATIONAL PURPOSES ONLY AND SHALL NOT BE USED FOR ANY PURPOSE OTHER THAN THAT AGREED UPON BETWEEN TOSHIBA AND THE CUSTOMER. THESE SITE PLANS ARE NOT TO BE USED FOR CONSTRUCTION PURPOSES.											
DATE:	06-19-13										
SCALE:	NOT TO SCALE										
PLANNER:	C.B.S.										
QUOTE:	N/A										
PROJECT NO.	130013698CTF										
E4											

FOR REFERENCE ONLY. NOT TO BE USED FOR CONSTRUCTION PURPOSES.